



**Aqualife
Products®**

**EXCELLENCE IN WATER TREATMENT
& LOW COST MAINTENANCE SYSTEMS**



**Quality Assurance
Guaranteed**

New AquaClave MARK 3 RO P5

It just keeps getting better!

- **Leak Detector - easy to use, hassle free and EFFECTIVE**
- **USA Manufactured housings - arguably world's best SECURITY**
- **UPGRADED TRANSFORMER - electrics Australian designed**

The NEW AquaClave Mark 3 remains a multi-stage system relying on five processes to deliver water to steam sterilisation equipment. The new model includes the automatic flushing system PLUS a leak detector, American housings and an improved power supply!

Meticulous attention to all stages of processes of the water includes 1 micron sediment removal and selection of the granular activated carbon specification.



Reverse Osmosis - The reverse osmosis *low production* AquaClave membrane has been engineered to produce exceptionally high quality in the most challenging water quality environment. To monitor performance a Total Dissolved Solids meter is supplied as part of the product.

The two stage polishing system will remove remaining organic and non-organic residues thereby finishing the water to meet all autoclave operating standards.

Commercial Advantages - The availability of this high quality water ON TAP avoids the purchasing of bottled distilled water, avoids storage challenges, minimises autoclave maintenance costs and reduces the damage to instruments from tap water rinsing makes buying an **AquaClave Mark 3** an intelligent decision.

Please note that this system is only suitable for use on disinfected potable water supplies from your local water authority. For areas where tank or other water is required please consult Aqualife. To achieve the water quality indicated above the system must be operated in accordance with the operating and maintenance instructions. On site verification of the water quality is required for municipal water supplies.

Dimensions: 46cm W x 54cm H x 23cm D